



1 12 December 2011
2 EMA/CVMP/VICH/463/2002-Draft 2
3 Committee for Medicinal Products for Veterinary Use (CVMP)

4 **VICH GL34: Biologicals: testing for the detection of**
5 **mycoplasma contamination**
6 Draft

Draft agreed by IWP	May 2002
Adoption by CVMP for release for consultation	15 May 2002
End of consultation (deadline for comments)	30 April 2003
Adoption by CVMP for release for further ¹ consultation	9 December 2010
End of consultation (deadline for comments)	31 March 2011
Draft 2 agreed by VICH Steering Committee	17 November 2011
Adoption by CVMP for release for consultation	8 December 2011
End of consultation (deadline for comments)	12 March 2012

7
8
9

Comments should be provided using this [template](#). The completed comments form should be sent to vet-guidelines@ema.europa.eu

10

¹ Draft VICH GL34 was published for consultation first in 2002. Following a 12 months it was agreed to suspend the consultation to wait for the testing of reference strains. The consultation was re-opened to allow submission of any further comments.





VICH GL34 (BIOLOGICALS: MYCOPLASMA)

November 2011

For consultation at Step 4 - Draft 2

11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26

TESTING FOR THE DETECTION OF
MYCOPLASMA CONTAMINATION

27
28
29
30
31
32
33
34
35
36
37
38
39
40
41

Recommended for Consultation
at Step 4 of the VICH Process
on 17 November 2011
by the VICH Steering Committee

THIS GUIDELINE HAS BEEN DEVELOPED BY THE APPROPRIATE VICH EXPERT WORKING GROUP AND IS SUBJECT TO CONSULTATION BY THE PARTIES, IN ACCORDANCE WITH THE VICH PROCESS. AT STEP 7 OF THE PROCESS THE FINAL DRAFT WILL BE RECOMMENDED FOR ADOPTION TO THE REGULATORY BODIES OF THE EUROPEAN UNION, JAPAN AND THE USA.

46 1. INTRODUCTION

47 1.1. Objective of the guideline

48 This VICH (International Cooperation on Harmonization of Technical Requirements for
49 Registration of Veterinary Medicinal Products) guideline is intended to facilitate the
50 harmonized licensing of new products for veterinary use. It is important that biological
51 products for veterinary use are free of contamination with Mycoplasmas to help assure
52 consistency of production and final product safety. Mycoplasma contaminants may be
53 introduced into cell culture and *in ovo* origin biological products through the master
54 seeds, the master cell seed (stock), starting materials of animal origin, and in
55 processing of biological materials during passage and product assembly. Therefore it
56 is necessary to demonstrate through testing that Mycoplasmas are not present, within
57 the limits of the test, in the final product, working seeds and cells and harvests, and
58 starting materials such as the master seed, master cell seed, and ingredients of animal
59 origin. This guideline establishes stages of manufacture to be tested and test
60 procedures to detect the presence of Mycoplasma contamination. It will provide a
61 unified standard that will facilitate the mutual acceptance of test data by the relevant
62 regulatory authorities. Methods proven equivalent to the guideline method by
63 scientifically accepted criteria could also be acceptable.

64 1.2. Background

65 The present methods for testing for Mycoplasma contamination are described in the
66 Japanese “Minimum requirements of biological products for animal use (2002)”, the
67 European Pharmacopoeia (7th Edition, 2011, 2.6.7), and the United States Code of
68 Federal Regulations, Title 9, 113.28. These requirements are all similar in that they
69 require testing for Mycoplasma contamination using a broth and agar technique. The
70 requirements do however differ in the specifics of these broth and agar tests as well as
71 other alternative test methods that are required or approved for use in detecting
72 Mycoplasma contamination.

73 1.3. Scope of guideline

74 This guideline describes the manner in which tests conducted to detect the presence
75 of Mycoplasma contamination in cell culture and *in ovo* origin biological products for
76 veterinary use shall be done to assure the absence of Mycoplasma contamination.
77 Tests on master seeds, master cell seeds (stocks), working seeds and cells,
78 ingredients of animal origin, harvests and live final vaccine and harvests for killed
79 products are included. Bacterial products which grow in the mycoplasma test media
80 and products for which mycoplasma contamination risk has been addressed through a
81 validated mycoplasma inactivation procedure will be considered outside this guideline.
82 The absence of mycoplasma contamination in eggs used for production is controlled
83 by appropriate testing of the flock, which is not covered by this guideline.

84 1.4. Test Methods

85 The guideline describes two test methods: 1) expansion in broth culture and detection
86 by colony formation on nutrient agar plates; and 2) expansion in cell culture and
87 characteristic fluorescent staining of DNA (a technique capable of detecting non-
88 cultivatable strains).

89
90 A third test methodology, nucleic acid amplification (NAT), is acknowledged, but is not
91 included in this guideline. The use of validated NAT techniques is currently approved
92 or under consideration by regulatory authorities for more rapid detection confirmation,

93 and strain identification. Appropriately validated NAT techniques may be used as an
 94 alternative to the broth/agar culture method and/or the indicator cell culture method
 95 provided the NAT test is shown to be at least equivalent in detection limit to the test
 96 methods in this guideline. A sample that tests positive by NAT may be directly
 97 considered unsuitable for use. If confirmation of the presence of living mycoplasma in
 98 the material under test is needed the broth/agar culture method or the indicator cell
 99 culture method should be used. Evaluation of NAT method use is encouraged in
 100 parallel testing to further develop, compare, and refine the technique for possible
 101 inclusion in future versions of this guideline.

102 2. GUIDELINE FOR TESTING FOR MYCOPLASMA 103 CONTAMINATION

104 2.1. General test procedures for detecting Mycoplasma contamination

105 The culture method using broth and agar is the fundamental method of Mycoplasma
 106 detection. A solid and liquid media culture method shall be used to test harvests or
 107 final batches of vaccine, and ingredients of animal origin. Master seed, master cell
 108 seed (stock), and working seed and cell lots shall be tested using both a solid and
 109 liquid media culture method and an indicator cell culture method with DNA stain.
 110 Should either method result in a positive test for mycoplasma the sample is considered
 111 positive and is unsuitable for use.
 112

Material	Broth & Agar Culture	DNA Stain
Master Seed & Master Cell Seed	Required	Required
Working Seed & Working Cell Seed	Required	Required
Ingredient of Animal Origin ^{1,2}	Required	
Harvest	When testing required ³	
Final Product	When testing required ³	

113 ¹ Excluding eggs

114 ² Unless a validated mycoplasma inactivation procedure has been applied

115 ³ The competent authorities require testing of different combinations of harvests and final
 116 product.

117 2.2. Culture test system validation

118 The culture method should be carried out to validate the detection limit of a
 119 laboratory's mycoplasma detection method. A sufficient number of both solid and
 120 liquid media shall be used to insure the growth of a low level of the following 5 strains
 121 of mycoplasmas.

122 *Acholeplasma laidlawii*

123 *Mycoplasma hyorhinitis*

124 *Mycoplasma orale*

125 *Mycoplasma synoviae*

126 *Mycoplasma fermentans*

127

128 The species were selected to reflect a range (within a practical number) of antibiotic
 129 sensitivity (to detect inhibition of mycoplasma growth in the assay), fastidiousness,
 130 rapidity of growth, likelihood of being a contaminant, and pathogenicity in avian or
 131 mammalian target species. *Acholeplasma laidlawii* is a common cell culture
 132 contaminant of animal and possibly environmental origin. *Mycoplasma hyorhinitis* is
 133 fastidious, a common cell culture contaminant of animal origin, and a mammalian
 134 pathogen. *Mycoplasma orale* is antibiotic sensitive and is a common cell culture

135 contaminant of human origin. *Mycoplasma synoviae* is fastidious (having a
 136 nicotinamide-adenine-dinucleotide [DPN, NAD] and cysteine requirement) and is an
 137 avian pathogen. *Mycoplasma fermentans* is a slow-growing organism and a common
 138 cell culture contaminant of human origin.

139
 140 References of the strains used to validate the laboratory mycoplasma contamination
 141 culture test system should be of low passage level (15 or less), and identified relative
 142 to type culture isolates, (see Appendix 3.2 for further information on reference strains).
 143 The reference strains used to validate the culture test system will be appropriate to the
 144 products tested (see table). Validation for *M. synoviae* is required when materials of
 145 avian origin are used at any stage in development and production. Validation for *M.*
 146 *hyorhinis* and *A. laidlawii* is required when materials of mammalian origin are used at
 147 any stage in development and production. Validation for *M. orale* is required when an
 148 antibiotic has been used at any stage in development and production. Reference
 149 Preparations shall be used to validate each production lot of broth and agar. At least
 150 one reference strain must be used as a control with each test.

151
 152 **Required Reference Organisms by: product type; test method, and presence of**
 153 **antibiotics**

154

Vaccine type Antibiotic content Test Method	<i>A.</i> <i>laidlawii</i>	<i>M.</i> <i>orale</i>	<i>M.</i> <i>hyorhinis</i>	<i>M.</i> <i>synoviae</i>	<i>M.</i> <i>fermentans</i>
Avian <i>in ovo</i> origin vaccine Without Antibiotics Broth/Agar Method				X	X
Avian <i>in ovo</i> origin vaccine With Antibiotics Broth/Agar Method		X		X	X
Avian cell culture origin vaccine Without Antibiotics Broth/Agar Method	X			X	X
Avian cell culture origin vaccine With Antibiotics Broth/Agar Method	X	X		X	X
Mammalian cell culture origin vaccine Without Antibiotics Broth/Agar Method	X		X		X
Mammalian cell culture origin vaccine With Antibiotics Broth/Agar Method	X	X	X		X
Vaccine Without Antibiotics DNA Staining Method		X	X		
Vaccine With Antibiotics DNA Staining Method		X	X		

155

156 **2.3. Culture Method**

157 **2.3.1. Incubation conditions**

158 Incubate the broth culture medium or media in tightly stoppered containers in
159 air. Incubate all agar plates under microaerophilic conditions (nitrogen
160 containing 5-10% CO₂). For the solid medium or media, maintain an
161 atmosphere of adequate humidity to prevent desiccation of the agar surface.

162 **2.3.2. Nutritive properties of a new batch of medium**

163 Each new lot (batch) of medium must be tested for the nutritive properties using
164 references specified above in **Section 2.2**. Each testing laboratory must
165 determine the inoculum for each of their references that will contain a low level
166 (not more than 100 CFU). Inoculate the solid medium with a low level (not
167 more than 100 CFU) per 60 mm plate and per 100 ml container of broth
168 medium. Use at least one agar plate and broth container for each reference.
169 Incubate the agar and broth media and make subcultures from the broth onto
170 agar at the specified intervals. The agar medium batch complies with the test
171 for nutritive properties if for all the references specified, growth obtained does
172 not differ by a factor greater than 5 from the value calculated with respect to the
173 inoculum. The broth complies if Mycoplasma growth on those agar plates
174 subcultured from the broth is achieved for each reference specified. Media
175 formulations found effective are included in **Appendix 3.1** of this guideline.

176 **2.3.3. Inhibitory substances**

177 Carry out the test for nutritive properties in the presence and absence of the
178 material to be tested at the time of prelicense and whenever there is a change
179 in the production method that may affect the detection of mycoplasmas. If
180 growth of the references occurs more than one sub-culture sooner without the
181 test material than with the test material, or if plates directly inoculated with the
182 test material have less than one-fifth of the colonies of those directly inoculated
183 without the test material, the test material contains inhibitory substances.
184 These substances must be neutralized or their effect otherwise countered, e.g.,
185 through passage in substrates not containing inhibitors or dilution in a larger
186 volume of medium, before the test for mycoplasma contamination is carried
187 out. For the dilution technique, larger medium volumes may be used or the
188 inoculum volume may be divided among multiple 100 ml flasks. The
189 effectiveness of the neutralization or other process is confirmed by repeating
190 the test for inhibitory substances after neutralization.

191 **2.3.4. Test method**

192 **2.3.4.1.** The amount of inoculum for each plate of solid medium is 0.2 ml of
193 product to be examined. When an assay for mycoplasma concerns
194 master and working seeds, master and working cells, and ingredients
195 of animal origin a volume of not less than 10 ml of undiluted sample
196 shall be tested in each liquid medium. The volume of final product to
197 be tested in each liquid medium shall be as required by the regulatory
198 authority issuing the marketing authorization. These are currently not
199 less than 1 ml in Japan and the US and not less than 10 ml in the EU.
200 Incubate the agar plates at 35°C to 38°C, microaerophilically, for 10-14
201 days in an atmosphere of adequate humidity to prevent desiccation of
202 the surface. Incubate the liquid media at 35° C to 38° C in tightly
203 stoppered containers in air for 20-21 days. At the same time incubate

204 an uninoculated 100 ml portion of each liquid medium and agar plates
205 as a negative control. If any significant pH change occurs upon the
206 addition of the product to be examined (this should be determined at
207 the time of prelicense), the liquid medium shall be restored to its
208 original pH value by the addition of a solution of either sodium
209 hydroxide or hydrochloric acid. Between the 2nd and 4th day after
210 inoculation, subculture each liquid culture by inoculating at least 1
211 plate of each solid medium with 0.2 ml and incubate them at 35°C to
212 38°C microaerophilically for 10-14 days. Repeat the procedure
213 between the 6th and 8th day, again between the 13th and 15th day, and
214 again between the 19th and 21st day of the test. Incubate those agar
215 plates inoculated on day 19, 20, or 21 for 7 days. Observe the liquid
216 medium or media every 2 or 3 days and if a color change occurs,
217 subculture. Color change detection requires the addition of phenol red
218 to the media.
219

220 **2.3.4.2.** If the liquid medium or media shows bacterial or fungal contamination,
221 repeat the test. If it is not possible to read at least one plate per
222 inoculation day, the test must be repeated.
223

224 **2.3.4.3.** Include in the test, positive controls prepared by inoculating a low level
225 (not more than 100 CFU) of at least one of the reference species onto
226 the agar plates and into the broth medium or media. If the test is run
227 on a routine basis, the control species should be rotated on a regular
228 basis. This control shall be used in each test conducted with a
229 medium that has been validated for nutritive properties using
230 references determined by the types of products being tested as
231 specified in **Section 2.2** of this guideline.
232

233 **2.3.5. Judgment of the culture method**

234 At the end of the incubation period, examine all the inoculated solid media
235 microscopically for the presence of mycoplasma colonies. The product is
236 negative for Mycoplasma contamination if the growth of typical Mycoplasma
237 colonies has not occurred on any of the inoculated solid media. If growth of
238 typical Mycoplasma colonies has occurred on any of the solid media, the test
239 and sample tested are considered positive for Mycoplasma contamination. The
240 test is invalid if the positive controls do not show growth of mycoplasma on at
241 least one subculture plate or the negative controls are positive for mycoplasma
242 contamination. If either of the controls is invalid the test must be repeated. If
243 suspect colonies are observed, confirmation of mycoplasma contamination may
244 be accomplished using an appropriate and validated method.

245 **2.4. Indicator cell culture method**

246 Cell cultures are stained with a fluorescent dye that binds to DNA. Mycoplasmas are
247 detected by their characteristic particulate or filamentous pattern of fluorescence on
248 the cell surface, and if contamination is heavy, in the surrounding areas. Mitochondria
249 in the cytoplasm may be stained, but may be differentiated from mycoplasma.

250 **2.4.1. Validation of the indicator cell culture method**

251 Using a VERO or other equivalent in efficiency indicator cell culture substrate,
252 validate the procedure using an inoculum not more than 100 CFU or CFU-like
253 micro-organisms of appropriate references of *M. hyorhinis* and *M. orale*. Both

254 references must be positive when stained with the DNA stain at the end of the
255 test.

256
257 If for viral, etc., suspensions the interpretation of results is affected by
258 cytopathic effects, the virus may be neutralized using a specific antiserum that
259 has no inhibitory effects on mycoplasmas, or an alternative cell culture
260 substrate that does not allow the growth of the virus may be used. To
261 demonstrate the absence of inhibitory effects of serum, carry out the positive
262 control tests in the presence of neutralizing antiserum. Antiserum lots may be
263 qualified once rather than at use.

264 **2.4.2. Test method**

265 **2.4.2.1.** Seed the indicator cell culture at a suitable density that will yield
266 confluence of the cells after 3 days of growth (example: 2×10^4 to $2 \times$
267 10^5 cells per ml, 4×10^3 to 2.5×10^4 cells/cm²) in a cell culture vessel
268 of not less than 25 cm². The indicator cell culture should be sub-
269 cultured without antibiotic prior to use. Inoculate 1 ml of the sample to
270 be examined into the cell culture vessel and incubate at 35° C to 38° C.

271
272 **2.4.2.2.** After at least 3 days of incubation and the cells have grown to
273 confluence, make a subculture onto cover slips in suitable containers
274 or on some other surface (chambered slides) suitable for the test
275 procedure. Seed the cells in the second subculture at a low density so
276 that they reach only 50% confluence after 3-5 days of incubation.
277 Complete confluence must be avoided because it impairs visualization
278 of mycoplasmas after staining.

279
280 **2.4.2.3.** Remove medium from cover slips or chambered slides. Rinse the
281 monolayer of indicator cells with phosphate buffered saline (PBS) and
282 then fix with glacial acetic acid/methanol (1 to 3) or some other
283 suitable fixing solution.

284
285 **2.4.2.4.** Remove the fixing solution and discard. Wash the fixing solution with
286 sterile water and dry slides completely if they are to be stained more
287 than one hour later.

288
289 **2.4.2.5.** Add a suitable fluorescent dye that binds to DNA such as bisbenzimidazole
290 stain (Hoechst compound 33258, bisbenzimidazole, 5 ug/L) and allow
291 to stain for a suitable time.

292
293 **2.4.2.6.** Remove the stain and rinse the monolayer with water. Mount the
294 cover slips if applicable and examine the slides by fluorescence (for
295 bisbenzimidazole stain use a 330 nm/380 nm excitation filter, LP 440 nm
296 barrier filter) at 400 X magnification or greater.

297
298 **2.4.2.7.** Compare the microscopic appearance of the test cultures with that of
299 the negative and reference controls, examining for extranuclear
300 fluorescence. Mycoplasmas produce pinpoints or filaments over the
301 indicator cell's cytoplasm. They may also produce pinpoints and
302 filaments in the intercellular spaces. Multiple microscopic fields as
303 validated should be examined.

304 **2.4.3. Judgment of the indicator cell culture method**

305 The product being examined is negative for Mycoplasma contamination if there
306 is no evidence of pinpoint or filaments of extranuclear fluorescence. If the
307 slides inoculated with the product contain evidence of pinpoint or extranuclear
308 fluorescence indicative of Mycoplasma the test and sample tested are
309 considered positive for Mycoplasma contamination. The test is invalid if the
310 positive controls do not show the presence of the appropriate extranuclear
311 fluorescence of the reference organisms or the negative cell controls contain
312 extranuclear fluorescence. If either of the controls is invalid the test must be
313 repeated.

314 3. APPENDICES

315

316 **3.1. Regional examples of suitable broth and agar formulations**

317

318 **9 CFR Mycoplasma Broth**

319	Heart Infusion Broth	62.5 g
320	Proteose Peptone #3	25.0 g
321	Yeast Extract	12.5 ml
322	1 % Thallium Acetate	62.5 ml
323	1 % Tetrazolium Chloride	13.75 ml
324	Penicillin (100,000 units/ml)	12.5 ml
325	Heat inactivated Horse Serum	250 ml
326	H ₂ O	2425 ml

327

328 Mix all ingredients well and adjust pH to 7.9 with 10 Normal NaOH.

329 Filter sterilize through a 0.2 µ filter. Dispense into sterile test vessels.

330 Add DPN/L-Cysteine solution before use, 2 ml/100 ml of broth.

331

332 **9 CFR Mycoplasma Agar**

333	Heart Infusion Agar	25 g
334	Heart Infusion Broth	10 g
335	Proteose Peptone #3	10 g
336	1% Thallium Acetate	25 ml
337	H ₂ O	995 ml
338	Heat Inactivated Horse Serum	126 ml
339	Yeast Extract	5 ml
340	Penicillin (100,000 units/ml)	5.2 ml
341	DPN/L-Cysteine	21 ml

342

343 Combine heart infusion agar, heart infusion broth, proteose peptone #3, Thallium acetate,
344 and H₂O.

345 Mix and bring to boil, then cool. Adjust the pH to 7.9 with 10 Normal NaOH.

346 Autoclave 20 min. at 121° C. Cool in water bath to 56° C.

347 Aseptically add: horse serum, yeast extract, Penicillin, and DPN/L-Cysteine.

348 Dispense 12 ml into each 15 X 60 mm petri dish.

349

350 **DPN/L-Cysteine solution**

351	Nicotiamide-adenine-dinucleotide (DPN, NAD)	5 g
352	Q.S. with H ₂ O to	500 ml

353

354	L-Cysteine	5 g
355	Q.S. with H ₂ O to	500 ml

356

357 Mix each chemical separately until dissolved.

358 Mix the two solutions and filter sterilize.

359

360 **Japanese Liquid Medium for Mycoplasma**

361	Basal Medium	
362	50 % w/v Bovine Cardiac Muscle Extract	100 ml
363	Meat Peptone	10 g

364	Sodium Chloride	5 g
365	Glucose	1 g
366	Sodium L-glutamate	0.1 g
367	L-arginine hydrochloride	1 g
368	H ₂ O	QS to 1000 ml
369		
370	Filter sterilize through 0.22 μ membrane filter or sterilize at 121° C for 15 min.	
371	Adjust the pH of the medium after sterilization to 7.2-7.4.	
372		
373	Additives for 77 ml of the Basal medium;	
374	Horse Serum	10 ml
375	Inactivated Porcine Serum	5 ml
376	25 % w/v Fresh Yeast Extract	5 ml
377	1 % w/v β-NAD (oxidized)	1 ml
378	1 % w/v L-cysteine HCL (1 H ₂ O)	1 ml
379	0.2 % w/v phenol red	1 ml
380		
381	Previously filter sterilize the additives and aseptically add to the sterilized basal	
382	medium. The additives which can be sterilized by high pressure can be autoclaved.	
383	Penicillin G potassium, 500 units/ml of the medium, or Thallium acetate, 0.02 % w/v,	
384	can be added.	
385		
386	Japanese Agar Medium for Mycoplasma	
387	Basal Medium	78 ml
388	Agar	1 g
389		
390	Sterilize by autoclaving 121°C for 15 min.	
391		
392	Additives:	
393	Horse Serum	10 ml
394	Inactivated Porcine Serum	5 ml
395	25 % w/v fresh yeast extract	5 ml
396	1 % w/v β-NAD (oxidized)	1 ml
397	1 % w/v L-cysteine HCl (1 H ₂ O)	1 ml
398		
399	Penicillin G potassium, 500 units per ml of medium, or thallium acetate, 0.02 % w/v can be	
400	added.	
401	Add the additives to basal/agar medium which has been liquefied by heating, and divide into	
402	sterile petri dishes, 45-55 mm. Cool and allow to solidify.	
403		
404	EP Hayflick media (Recommended media for the general detection of <i>mycoplasmas</i>)	
405	Liquid Medium:	
406	Beef Heart Infusion Broth (1)	90 ml
407	Horse Serum (unheated)	20 ml
408	Yeast Extract (250 g/L)	10 ml
409	Phenol Red (0.6 g/L solution)	5 ml
410	Penicillin (20,000 I.U. per ml)	0.25 ml
411	Deoxyribonucleic acid (2 g/L solution)	1.2 ml
412		
413	Adjust to pH 7.8.	
414		
415	Solid Medium:	

416 Prepare as described for the liquid medium above but replace beef heart infusion broth
417 with beef heart infusion agar containing 15 g/L of agar.

418

419 **EP Frey media (Recommended Media for the detection of *M. synoviae*)**

420 Liquid Medium:

421	Beef Heart Infusion Broth (1)	90 ml
422	Essential Vitamins (2)	0.025 ml
423	Glucose monohydrate (500 g/L solution)	2 ml
424	Swine serum (inactivated at 56°C for 30 min.)	12 ml
425	β-Nicotinamide adenine dinucleotide (10 g/L solution)	1 ml
426	Cysteine hydrochloride (10 g/L solution)	1 ml
427	Phenol Red (0.6 g/L solution)	5 ml
428	Penicillin (20,000 I.U. per ml)	0.25 ml

429

430 Mix the solutions of β-nicotinamide adenine dinucleotide and cysteine hydrochloride
431 and after 10 minutes, add the other ingredients. Adjust pH to 7.8.

432

433 Solid Medium:

434	Beef Heart Infusion Broth (1)	90 ml
435	Ionagar (3)	1.4 g

436

437 Adjust pH to 7.8, and sterilize by autoclaving, then add:

438	Essential Vitamins (2)	0.025 ml
439	Glucose monohydrate (500 g/L solution)	2 ml
440	Swine serum (unheated)	12 ml
441	β-Nicotinamide adenine dinucleotide (10g/L solution)	1 ml
442	Cysteine hydrochloride (10 g/L solutions)	1 ml
443	Phenol Red (0.6 g/L solution)	5 ml
444	Penicillin (20,000 I.U. per ml)	0.25 ml

445

446 **EP Friis media (Recommended Media for the Detection of Non-avian Mycoplasmas)**

447 Liquid Medium:

448	Hank's Balanced Salt Solution (modified) (4)	800 ml
449	H ₂ O	67 ml
450	Brain Heart Infusion (5)	135 ml
451	PPLO Broth	248 ml
452	Yeast Extract (170 g/L)	60 ml
453	Bacitracin	250 mg
454	Meticillin	250 mg
455	Phenol Red (5 g/L)	4.5 ml
456	Horse Serum	165 ml
457	Swine Serum	165 ml

458

459 Adjust the pH to 7.40-7.45

460

461 Solid Medium:

462	Hank's Balanced Salt Solution (modified) (4)	200 ml
463	DEAE-dextran	200 mg
464	Ionagar (3)	15.65 g

465

466 Mix well and sterilize by autoclaving. Cool to 100° C. Add this to 1740 ml of the liquid
467 medium described above.

468

469 **EP Media Sub parts**

470 (1) Beef Heart Infusion Broth

471	Beef Heart (for preparation of the infusion)	500 g
472	Peptone	10 g
473	Sodium Chloride	5 g
474	H ₂ O	QS to 1000 ml

475
476 Sterilize by autoclaving.

477

478 (2) Essential Vitamins

479	Biotin	100 mg
480	Calcium pantothenate	100 mg
481	Choline chloride	100 mg
482	Folic acid	100 mg
483	<i>l</i> -Inositol	200 mg
484	Nicotinamide	100 mg
485	Pyridoxal hydrochloride	100 mg
486	Riboflavin	10 mg
487	Thiamine hydrochloride	100 mg
488	H ₂ O	QS to 1000 ml

489

490 (3) Ionagar

491 A highly refined agar for use in microbiology and immunology, prepared by an ion-
492 exchange procedure which results in a product having superior purity, clarity, and gel
493 strength.

494 It contains approximately:

495	H ₂ O	12.2 %
496		
497	Ash	1.5 %
498	Acid insoluble ash	0.2 %
499	Chlorine	0.0 %
500	Phosphate (calculated as P ₂ O ₅)	0.3 %
501	Total Nitrogen	0.3 %
502	Copper	8 ppm
503	Iron	170 ppm
504	Calcium	0.28 %
505	Magnesium	0.32 %

506

507 (4) Hank's Balanced Salt Solution (modified)

508	Sodium chloride	6.4 g
509	Potassium chloride	0.32 g
510	Magnesium sulphate heptahydrate	0.08 g
511	Magnesium chloride hexahydrate	0.08 g
512	Calcium chloride, anhydrous	0.112 g
513	Disodium hydrogen phosphate dihydrate	0.0596 g
514	Potassium dihydrogen phosphate, anhydrous	0.048 g
515	H ₂ O	QS to 800 ml

516

517 (5) Brain heart infusion

518	Calf brain infusion	200 g
519	Beef heart infusion	250 g
520	Proteose peptone	10 g
521	Glucose	2 g

522	Sodium chloride	5 g
523	Disodium hydrogen phosphate, anhydrous	2.5 g
524	H ₂ O	QS to 1000 ml
525		
526	(6) PPLO broth	
527	Beef heart infusion	50 g
528	Peptone	10 g
529	Sodium chloride	5 g
530	H ₂ O	QS to 1000 ml
531		

532 **Bisbenzimidazole stain solution for DNA Staining**

533 Hoechst compound 33258 (bisbenzimidazole), 5 µg per liter of buffered aqueous
534 solution.

535
536 Note: The solution should be protected from light.

537 **3.2. Mycoplasma References**

538 Standardization of testing between laboratories and between regions would be
539 enhanced by use of references common within or between regions, This has been
540 shown to be presently impractical due to the difficulty of producing consistent batches
541 of lyophilized references and the shipping issues with frozen references. Therefore,
542 regions or laboratories may use their own references, providing that they are of low
543 passage level (15 or less), identified relative to type culture isolates, stable, and
544 appropriately validated as suitable for use in the context of this guideline. It is strongly
545 recommended to include in the validation of detection limit a comparison to the
546 EDQM reference strains (described below) for international recognition. Regions or
547 laboratories may produce their own validated references, or may acquire commonly
548 available and appropriately validated references, such as the following produced by
549 the EDQM.

550
551 The 5 strains of Mycoplasma listed in **Section 2.3** were isolated by laboratories of the
552 European Union and donated to the European Directorate for the Quality of Medicines
553 and HealthCare (EDQM). EDQM produced a sufficient quantity of these frozen
554 references, and performed an intra-region EU validation/stability study (C. Milne, A.
555 Daas. Establishment of European Pharmacopoeia Mycoplasma Reference Strains.
556 *Pharmeuropa Bio* 2006(1):57-72). Completion of further validation studies by the
557 regulatory agencies and the industries in Japan, USA, and Canada confirmed that the
558 strains are very suitable for use in the context of this guideline. (VICH Collaborative
559 Study on the Ph. Eur. Mycoplasma Reference Strains: EDQM Report Compiling and
560 Analysing the Data Set for the VICH Collaborative Study on the European
561 Pharmacopoeia Mycoplasma Reference Strains, EDQM Administrator Representative,
562 C. Milne, 2010.) The BQMEWG commends the Staff at the EDQM for their efforts and
563 perseverance in producing and validating these very excellent references.

564
565
566 For DNA staining validation, the following strains may also prove useful:

567 *M. hyorhinis* -- ATCC 29052

568 *M. orale* -- ATCC 23714

569

570 **3.3. Glossary**

571
572
573
574
575
576
577
578
579
580
581
582
583
584
585
586
587
588
589
590
591
592
593
594
595
596
597
598
599
600
601
602
603
604
605
606
607
608
609
610
611
612
613
614
615
616
617
618
619
620
621
622
623

Batch (lot, serial) of starting material of animal origin

The total quantity of homogenous material (e.g., cells, serum) identified by a unique serial number.

Cell-seed system

A system whereby successive final lots (batches) of a product are manufactured by culture in cells derived from the same master cell seed. A number of containers from the master cell seed are used to prepare a working cell seed.

Cell lines

Cultures of cells >10 passages or subcultures from the tissue of origin and having a high capacity for multiplication *in-vitro*.

Final product, batch, lot, or serial

A collection of closed, final containers or other final dosage units that are expected to be homogeneous and equivalent with respect to risk of contamination during filling or preparation of the final product. The dosage units are filled, or otherwise prepared, from the same final bulk product, freeze-dried together (if applicable) and closed in one continuous working session. They bear a distinctive number or code identifying the final lot (batch, serial). Where a final bulk product is filled and/or freeze-dried on several separate sessions, there results a related set of final lots (batches, serials) that are usually identified by the use of a common part in the distinctive number or code; these related final lots (batches, serials) are sometimes referred to as sub-batches, subserials, sub-lots or filling lots. For the purposes of mycoplasma testing, a single sub-batch may be considered representative of the batch.

Harvests

Material derived on one or more occasions from a single production culture inoculated with the same working seed lot (single harvest) or pooled material containing a single strain or type of micro-organism or antigen and derived from a number of eggs, cell culture containers, etc. that are processed at the same time (monovalent pooled harvest).

Master cell seed (stock)

A collection of aliquots of cells (primary or cell line) of a single passage level for use in the preparation of the product, distributed into containers in a single operation, processed together and stored in such a manner as to ensure uniformity and stability and to prevent contamination. Master cell seed is usually stored at temperatures of -70°C or lower.

Master seed

A collection of closed containers of a culture of micro-organisms of a single passage level used for the production of all batches of a designated veterinary biological product, distributed from a single bulk into containers and processed together in a single operation in such a manner as to ensure uniformity and stability and to prevent contamination.

Microaerophilic condition

A nitrogen atmosphere containing 5-10% carbon dioxide and sufficient humidity to prevent drying of the agar plates.

624 **Passage**
625 One transfer of cells or microorganisms followed by the normally used incubation
626 period for the cell or microorganism concerned.
627

628 **Primary cell cultures**
629 Primary cell cultures are cultures of cells essentially unchanged from those in the
630 animal tissues from which they have been prepared and being no more than 10 *in-vitro*
631 passages to the test level from the initial preparation from the animal tissue. The first
632 *in-vitro* cultivation is regarded as the first passage of the cells
633

634 **Seed-lot system:**
635 A system in which successive batches of a product are derived from the same master
636 seed virus. For routine production, a working seed virus may be prepared from the
637 master seed virus.
638

639 **Working cell seed (stock)**
640 A collection of aliquots of cells derived from the master cell seed and at the passage
641 level used in the preparation of production cell cultures. The working cell seed is
642 distributed into containers, processed and stored as described for master cell seed.
643 The term includes production cell seed.
644

645 **Working References**
646 A passage of the Reference strains of Mycoplasma produced in the testing laboratory
647 for use as controls to satisfy the reference requirements specified in this document.
648

649 **Working seed**
650 A collection of aliquots of a microorganism derived from the master seed virus and at
651 the passage level used in the preparation of product. Working seed virus is distributed
652 into containers and stored as described for master seed virus. The term includes
653 production seed.
654
655