

## Assessment of Drug Sensitivity in Hematopoietic Stem and Progenitor Cells From Acute Myelogenous Leukemia and Myelodysplastic Syndrome Ex Vivo

KATHERINE L.B. KNORR,<sup>a</sup> LAURA E. FINN,<sup>b</sup> B. DOUGLAS SMITH,<sup>c</sup> ALLAN D. HESS,<sup>c</sup> JAMES M. FORAN,<sup>b</sup> JUDITH E. KARP,<sup>c</sup> SCOTT H. KAUFMANN<sup>a,d</sup>

**Key Words.** Hematopoietic stem cells • Myeloid progenitor cells • Myelodysplastic syndrome • Acute myelogenous leukemia • MLN4924

### ABSTRACT

Current understanding suggests that malignant stem and progenitor cells must be reduced or eliminated for prolonged remissions in myeloid neoplasms such as acute myelogenous leukemia (AML) or myelodysplastic syndrome (MDS). Multicolor flow cytometry has been widely used to distinguish stem and myeloid progenitor cells from other populations in normal and malignant bone marrow. In this study, we present a method for assessing drug sensitivity in MDS and AML patient hematopoietic stem and myeloid progenitor cell populations ex vivo using the investigational Nedd8-activating enzyme inhibitor MLN4924 and standard-of-care agent cytarabine as examples. Utilizing a multicolor flow cytometry antibody panel for identification of hematopoietic stem cells, multipotent progenitors, common myeloid progenitors, granulocyte-monocyte progenitors, and megakaryocyte-erythroid progenitors present in mononuclear cell fractions isolated from bone marrow aspirates, we compare stem and progenitor cell counts after treatment for 24 hours with drug versus diluent. We demonstrate that MLN4924 exerts a cytotoxic effect on MDS and AML stem and progenitor cell populations, whereas cytarabine has more limited effects. Further application of this method for evaluating drug effects on these populations ex vivo and in vivo may inform rational design and selection of therapies in the clinical setting. *STEM CELLS TRANSLATIONAL MEDICINE* 2016;5:1–11

### SIGNIFICANCE

Despite the fact that malignant hematopoietic stem and progenitor cells play critical roles in the initiation, maintenance, and progression of acute myelogenous leukemia (AML) and myelodysplastic syndrome (MDS), information about the effects of investigational chemotherapeutic agents on these cells is limited. Using the Nedd8-activating enzyme inhibitor MLN4924 as an example, the authors describe a method for examining the impact of novel chemotherapeutics on these populations ex vivo, demonstrate that MLN4924 is generally more active against AML stem and progenitor cells than the standard-of-care agent cytarabine, and report that MLN4924 is also active against MDS stem and progenitor cells ex vivo.

### INTRODUCTION

Multicolor flow cytometry has been previously used to distinguish various stem, progenitor, and mature cell populations within the bone marrow [1, 2]. To our knowledge, this approach has not previously been applied to drug sensitivity testing either ex vivo or in vivo. As a result, information about the impact of either investigational agents or standard-of-care drugs on stem and progenitor cell subpopulations is limited.

At the same time, there is substantial evidence that malignant hematopoietic stem and progenitor cells contribute to disease initiation, maintenance, and progression, particularly in

acute myelogenous leukemia (AML) and myelodysplastic syndrome (MDS) [1–7]. These malignant hematopoietic stem and progenitor cell populations, which have characteristic surface antigen profiles, have been shown to produce disease upon engraftment, exhibit abnormal differentiation capacity, and outgrow nonmalignant cells [1, 3, 8–10]. Based on these observations, it is currently thought that targeting one or more of these malignant hematopoietic stem and progenitor populations might be necessary to achieve durable disease remissions in these disorders [2, 4, 5, 11–19]. Accordingly, an assay that assesses the impact of potential therapeutic agents on malignant versus normal

<sup>a</sup>Department of Molecular Pharmacology and Experimental Therapeutics and <sup>d</sup>Division of Oncology Research, Mayo Clinic, Rochester, Minnesota, USA;

<sup>b</sup>Division of Hematology and Oncology, Mayo Clinic, Jacksonville, Florida, USA;

<sup>c</sup>Division of Hematological Malignancies, Sidney Kimmel Cancer Center, Johns Hopkins University, Baltimore, Maryland, USA

Correspondence: Scott Kaufmann, M.D., Ph.D., Mayo Clinic, Gonda 19-121, 200 First Street Southwest, Rochester, Minnesota 55905, USA. Telephone: 507-284-8950; E-Mail: kaufmann.scott@mayo.edu

Received January 18, 2016; accepted for publication September 19, 2016.

©AlphaMed Press 1066-5099/2016/\$20.00/0

<http://dx.doi.org/10.5966/sctm.2016-0034>

hematopoietic stem and progenitor populations might help identify treatments that have unique activity in AML and/or MDS.

To address this issue, we have developed an *ex vivo* approach for assessment of drug sensitivity in hematopoietic stem and myeloid progenitor cell populations from AML and MDS patients. Using multicolor flow cytometry to assess malignant stem and progenitor cell populations, we compared the number of live hematopoietic stem cells (HSCs), multipotent progenitors (MPPs), common myeloid progenitors (CMPs), granulocyte-monocyte progenitors (GMPs), and megakaryocyte-erythroid progenitors (MEPs) in bone marrow mononuclear cell fractions cultured with drug versus diluent (control) for 24 hours *ex vivo*.

In this study, we describe the application of this approach using two different agents. Cytarabine, a nucleoside analog that is incorporated into DNA [20] and induces apoptosis [21], is widely used in AML therapy and is thought to be the single most active agent for the treatment of this group of disorders [22, 23]. Unfortunately, up to one third of AML patients fail to achieve a complete remission with cytarabine-based therapy, and many of those who do achieve a remission subsequently relapse [22, 23]. The second agent, MLN4924, is a first-in-class Nedd8-activating enzyme (NAE) inhibitor that prevents processing of the protein modifier Nedd8, which is required for activation of the Cullin-RING ligase (CRL) class of E3 ubiquitin ligases [24]. Without the NEDD8 modification, CRLs are inactive and their substrates accumulate, leading to cytotoxicity in neoplastic cells but largely sparing normal cells [25–27]. MLN4924 has shown promising clinical activity in myeloid malignancies, even in the relapsed or refractory setting [28].

Under the conditions of the *ex vivo* assay described in this study, MLN4924 is toxic to leukemic and myelodysplastic hematopoietic stem and progenitor cells, whereas cytarabine in general is not. The same assay shows that MLN4924 has limited toxicity in normal hematopoietic stem and progenitor cells, extending a recent report that assessed toxicity of MLN4924 in normal CD34<sup>+</sup> cells [29]. Thus, the strategy described here might help inform development of future AML and MDS therapies aimed at targeting malignant stem and progenitor cell populations while sparing normal hematopoietic progenitors.

## MATERIALS AND METHODS

### Chemicals and Reagents

Reagents were obtained from the following suppliers: MLN4924 (Active Biochem, Maplewood, NJ, <http://www.activebiochem.com>); enhanced chemiluminescence reagents and 16% paraformaldehyde (PFA; Thermo Fisher Scientific, Waltham, MA, <https://corporate.thermofisher.com>); cytarabine (AraC), heat-inactivated fetal bovine serum (FBS), dimethyl sulfoxide (DMSO), propidium iodide (PI), Triton X-100, and Histopaque 1077 (Sigma-Aldrich, St. Louis, MO, <https://www.sigmaaldrich.com>); and the broad spectrum caspase inhibitor N-(2-quinolyl)valyl-aspartyl-(2,6-difluorophenoxy)methyl ketone (Q-VD-OPH; SM Biochemicals, Anaheim, CA, <http://www.smbiochemicals.com>). Primary antibodies were purchased as follows: Cullin1, Bcl-x<sub>L</sub>, Mcl-1, Bim, and glyceraldehyde-3 phosphate dehydrogenase (Cell Signaling Technology, Beverly, MA, <https://www.cellsignal.com>); Bcl-2 (Dako, Carpinteria, CA, <http://www.dako.com>); Puma

(catalog no. sc-374223; Santa Cruz Biotechnology, Santa Cruz, CA, <https://www.scbt.com>); and Noxa (Enzo Life Sciences, Farmingdale, NY, <http://www.enzolifesciences.com>).

Antibodies and reagents for the flow cytometry panel were obtained as follows: APC-Annexin V, BV510-CD7, BV711-CD33, PE-Cy7-CD34, PE-CD47, FITC-CD123, Brilliant Violet Staining Buffer, Compensation Particles (Antimouse Igκ; BD Bioscience, Franklin Lakes, NJ, <http://www.bdbiosciences.com>); BV510-CD3, -CD10, -CD14, -CD15, -C19, -CD56, -CD64, AlexaFluor700-CD38, PacBlue-CD45, PerCp-Cy5.5-CD90, and Red Cell Lysis Buffer (BioLegend, San Diego, CA, <http://www.biolegend.com/>); APC-eFluor780-CD45RA (eBioscience, San Diego, CA, <http://www.ebioscience.com/>); and Aqua Fixable Live/Dead Stain (Thermo Fisher).

### Cells and Cell Culture

After institutional review board approval and patient consent, primary MDS and AML bone marrow aspirates were obtained as material in excess of that required for hematopathological examination during routine clinical care. Mononuclear cells from bone marrow aspirates were isolated using a standard Ficoll gradient procedure [30] and resuspended in RPMI 1640 medium containing 10% FBS (medium A) for subsequent experiments. Deidentified umbilical cord blood units containing normal stem and progenitor cells [31–34] collected from consenting donors and banked for research purposes were obtained from the Colorado Stem Cell Bank. Mononuclear cells were isolated from cord blood units using a previously published protocol [35] and resuspended in medium A.

### Detection of Apoptosis via PI Staining

Cells were exposed to increasing concentrations of MLN4924 in 0.1% DMSO for 24 hours, pelleted, and stained with PI in 0.1% (w/v) sodium citrate and 0.1% (w/v) Triton X-100 as previously described [36, 37]. After 30 minutes, samples were subjected to flow cytometry using a FACSCanto II (BD Bioscience) to assess PI fluorescence intensity. Analysis and quantification of the sub-G<sub>0</sub>/G<sub>1</sub> population were completed using FlowJo software (TreeStar, Ashland, OR, <http://www.flowjo.com>).


### Western Blotting

Cells exposed to diluent (0.1% DMSO) or increasing concentrations of MLN4924 for 24 hours were pelleted for protein isolation. Sample preparation, SDS-PAGE, transfer to nitrocellulose membrane, and subsequent protein detection by immunoblotting were performed as previously described [38]. These experiments were conducted in the presence of Q-VD-OPH, a broad spectrum caspase inhibitor [39], to limit cleavage or release of proteins and transcripts of interest (due to apoptosis and secondary to loss of membrane integrity) without affecting apoptotic events upstream of caspase activation.

### Quantification of Hematopoietic Stem and Progenitor Cells

Unsorted mononuclear cell fractions isolated from normal donor cord blood units or MDS and AML bone marrow aspirates were counted and plated at a density of  $2 \times 10^6$  cells in 3 ml medium A per well in 6-well plates.

Because the baseline number of stem and progenitor cells varies among cord blood samples and among bone marrow aspirates, each cord blood unit or bone marrow aspirate constitutes an independent experiment. Therefore, for each cord blood unit



**Figure 1.** Quantification of hematopoietic stem and progenitor cell populations in cord blood units and relative survival of populations in MLN4924-treated versus control samples. **(A):** Gating scheme for identification of hematopoietic stem and progenitor cell populations from initial CD34<sup>+</sup> CD45<sup>dim</sup> population. See text for details. **(B):** Antigen profiles for each population of interest. **(C):** Survival of CD34<sup>+</sup> CD45<sup>dim</sup> cord blood populations in MLN4924 or AraC compared with diluent control. Kruskal-Wallis test, 6 sample groups, *p* value = NS (.71). Error bars represent  $\pm$  SEM from four independent cord blood samples. Abbreviations: AraC, cytarabine; CMP, common myeloid progenitor; DMSO, dimethyl sulfoxide; GMP, granulocyte-monocyte progenitor; HSC, hematopoietic stem cell; MEP, megakaryocyte-erythroid progenitor; MPP, multipotent progenitor; NS, not significant.

or bone marrow aspirate, one well was treated with diluent (0.1% DMSO) to serve as the control sample for that experiment. Other wells were treated with various concentrations of MLN4924 (30–1,000 nM), consistent with previously published studies [25, 40]. Experiments with AML cells also included a well treated with 100 nM AraC, a concentration that reflects typical serum levels achieved with standard 7 + 3 chemotherapy regimens [41, 42].

After incubation for 24 hours with diluent or the indicated agent, samples were centrifuged at 150g for 5 minutes and washed once with cold PBS. Cells were then stained for 25 minutes at 4°C with Aqua Live/Dead stain and the fluorescent antibodies listed in the Chemicals and Reagents section. During this time, compensation beads were prepared in parallel according to manufacturer's instructions to establish appropriate flow cytometer settings before sample acquisition.

Labeled cells were washed twice in cold PBS, fixed in 2% PFA, and analyzed on an LSRII within 2 hours of preparation. In each experiment, approximately equal numbers of events were collected for each sample (diluent treated and drug treated), with  $0.7\text{--}1.0 \times 10^6$  events collected in various experiments. All analysis was performed using FlowJo software. Cell survival in the drug treated populations was expressed relative to the corresponding diluent treated population.

## RESULTS

### Method for Assessing Hematopoietic Stem and Progenitor Cell Survival in Ex Vivo Drug Sensitivity Testing

The present study was undertaken to compare the effects of cytarabine and the investigational agent MLN4924 in leukemic

stem and progenitor cells ex vivo. Normal cord blood mononuclear cells served two purposes in this study: (a) to ensure that the fluorochromes in the antibody panel could be compensated to establish a reliable gating scheme for quantification of stem and progenitor cells, and (b) to test whether MLN4924 or cytarabine was toxic to any of the normal populations of interest.

To eliminate dead and differentiated, lineage positive (Lin<sup>+</sup>) cells from analysis, a “dump gate” was established using Aqua Live/Dead stain (emission in the BV510 channel) and BV510 conjugated antibodies against antigens expressed on terminally differentiated cell types. The BV510<sup>-</sup> population consists of the live and undifferentiated, lineage negative (Lin<sup>-</sup>) cells that will be gated for analysis (Fig. 1A, far left panel). From the Live/Lin<sup>-</sup> gate, the CD34<sup>+</sup> CD45<sup>dim</sup> population containing all stem and progenitor cells was gated (Fig. 1A, arrow 1) and subsequently divided into stem and progenitor populations based on CD38 expression (Fig. 1A, arrow 2). CD34<sup>+</sup> CD38<sup>-</sup> cells (stem cells) were divided into HSCs and MPPs based on CD90 expression (Fig. 1A, arrow 3). CD34<sup>+</sup> CD38<sup>+</sup> cells (progenitors) were divided into CMPs, GMPs, and MEPs based on CD45RA and CD123 expression (Fig. 1A, arrow 4). The antigen profiles of each population of interest are summarized in Figure 1B. The relative survival of the bulk stem and progenitor cell population (CD34<sup>+</sup> CD45<sup>dim</sup>) in cord blood samples (*n* = 4) treated with MLN4924 or cytarabine did not differ significantly from the CD34<sup>+</sup> CD45<sup>dim</sup> population in the control samples exposed to diluent (Fig. 1C). These results are consistent with a recent report demonstrating lack of MLN4924 toxicity in normal human CD34<sup>+</sup> cells [29]. In the

**Table 1.** Cell counts in AML stem and progenitor populations

Patient	Sample	Total event count	Live/Lin <sup>-</sup> count	34 <sup>+</sup> 45 <sup>dim-a</sup>	Stem cells <sup>b</sup>	HSC	MPP	Myeloid progenitors <sup>c</sup>	CMF	GMP	MEP
1	DMSO	709,674	137,547	69,257	1,774	704	911	67,146	6,292	60,358	28
	30 nM MLN4924	705,327	136,210	68,247	2,002	840	975	65,874	6,604	58,750	30
	100 nM MLN4924	710,220	131,503	62,121	2,197	960	982	59,526	7,029	51,996	28
	300 nM MLN4924	708,009	119,324	53,481	2,244	853	1,130	50,872	6,265	44,140	33
	1,000 nM MLN4924	707,511	107,556	42,130	2,059	810	974	39,697	5,760	33,566	17
	DMSO	1,000,000	143,982	52,072	1,653	1,116	295	49,980	7,534	37,442	496
2	30 nM MLN4924	1,002,162	141,863	49,967	2,204	1,475	464	47,252	4,630	38,661	261
	100 nM MLN4924	1,000,000	132,661	42,000	2,060	1,404	439	39,456	3,676	32,229	220
	300 nM MLN4924	1,000,000	116,179	25,422	1,641	1,104	321	23,436	2,377	18,968	115
	1,000 nM MLN4924	1,006,825	100,017	17,469	1,481	963	292	16,107	1,828	12,703	72
	100 nM AraC	1,000,000	136,749	46,396	1,853	1,226	399	44,147	4,734	35,055	348
	DMSO	644,976	89,374	2,733	35	33	ND	2,672	2,509	43	62
3	1000 nM MLN4924	669,228	106,322	1,323	22	19	ND	1,284	1,184	47	22
	DMSO	378,757	29,619	10,651 (36% of Live/Lin <sup>-</sup> )	ND	ND	ND	ND	ND	ND	ND
	30 nM MLN4924	351,446	25,834	8,507 (33% of Live/Lin <sup>-</sup> )	878	288	19	19,523	4,183	13,484	288
	100 nM MLN4924	232,678	17,234	4,837 (28% of Live/Lin <sup>-</sup> )	720	239	25	16,440	3,867	11,067	269
	300 nM MLN4924	332,815	24,344	8,237 (34% of Live/Lin <sup>-</sup> )	466	128	16	12,626	3,052	8,353	245
	1,000 nM MLN4924	680,449	54,134	9,639	143	53	9	9,465	2,571	6,100	96
4	1,000 nM MLN4924	680,449	47,435	6,804	67	38	5	6,721	1,992	4,288	37
	100 nM AraC	680,449	39,292	14,803	742	262	15	13,949	3,350	9,218	262
	DMSO	1,006,383	184,582	88,150	27,505	1,983	1,166	58,551	3,946	54,075	8
	30 nM MLN4924	1,019,907	197,285	82,652	27,456	1,912	980	53,155	3,645	49,045	20
	100 nM MLN4924	1,007,811	172,364	69,633	22,244	1,392	762	45,748	2,793	42,660	5
	300 nM MLN4924	1,026,711	186,357	57,415	18,710	1,394	709	37,408	2,756	34,265	0
5	1,000 nM MLN4924	1,027,005	167,546	35,451	11,315	910	414	23,395	1,691	21,470	1
	100 nM AraC	1,014,468	177,369	75,359	28,389	2,999	1,166	45,102	4,353	39,932	3
	DMSO	1,000,000	150,817	35,963	9,610	4,852	3,133	25,469	3,126	21,889	73
	30 nM MLN4924	1,000,000	137,280	32,023	8,428	4,210	2,763	22,829	2,922	19,535	50
	100 nM MLN4924	1,000,000	122,138	27,201	7,347	3,573	2,429	19,231	2,664	16,212	72
	300 nM MLN4924	1,000,000	122,042	23,060	6,247	2,942	2,093	16,264	2,313	13,622	66
6	1,000 nM MLN4924	1,000,000	113,337	16,415	4,212	1,844	1,448	11,821	1,933	9,570	53
	100 nM AraC	1,000,000	119,818	25,864	7,573	4,003	2,335	17,606	2,282	14,936	73

Abbreviations: AML, acute myelogenous leukemia; AraC, cytarabine; CMF, common myeloid progenitor; DMSO, dimethyl sulfoxide; GMP, granulocyte-monocyte progenitor; HSC, hematopoietic stem cell; Lin<sup>-</sup>, lineage negative; MEP, megakaryocyte-erythroid progenitor; MPP, multipotent progenitor; ND, population not detected.

<sup>a</sup>CD34<sup>+</sup>CD45<sup>dim</sup>: starting population for analysis of all subpopulations (Fig. 1A).

<sup>b</sup>Stem cells: HSC + MPP.

<sup>c</sup>Myeloid progenitors: CMP + GMP + MEP.

present study, no decrease in survival was noted in normal HSCs, MPPs, CMPs, or GMPs during MLN4924 treatment (supplemental online Fig. 1). Although some toxicity was noted when normal MEPs were treated with MLN4924, it was less than the toxicity seen with cytarabine (supplemental online Fig. 1).

### MLN4924 Is Cytotoxic to De Novo and Secondary AML Hematopoietic Stem and Progenitor Cells

In contrast to the relative lack of effect of MLN4924 in normal stem and progenitor cells, MLN4924 preferentially diminished populations of leukemic stem and progenitor cells *ex vivo*. Using samples from seven AML patients (supplemental online Table 1), we compared surviving stem and progenitor cell populations in drug treated samples versus diluent control sample as illustrated in Figure 1. The number of events collected, cell count in the starting Live/Lin<sup>-</sup> gate used for analysis (i.e., the population used for analysis, as illustrated in the left panel of Fig. 1A), and the cell counts in each stem and progenitor population are reported in Table 1. Results of these analyses are summarized in Figure 2.

In *de novo* AML, MLN4924 treatment caused greater reduction of stem and progenitor cell populations than cytarabine did (Fig. 2, patients 1 through 5). However, two instances of resistance to MLN4924 were noted. First, decreased sensitivity to MLN4924 was noted in both stem (HSC and MPP) and progenitor (CMP and MEP) populations isolated from AML patient 1, who had FLT3 ITD<sup>+</sup> disease (Fig. 2, open triangles; Table 1). Only the GMP population in this AML sample was reduced by MLN4924. Second, in comparison with AraC sensitivity in the same clinical sample, AML patient 2 HSC and MPP populations appeared less sensitive to MLN4924 than the progenitor cell populations (Fig. 2, open squares). Specifically, reduction of HSC and MPP populations below the levels resulting from AraC treatment was not achieved until the MLN4924 concentration in this sample reached 300 nM. However, reduction in progenitor cell populations treated with only 30 nM MLN4924 was equal or greater to the reduction achieved with AraC.

### MLN4924 Is Also Cytotoxic to MDS and Post-MDS AML Hematopoietic Stem and Progenitor Cells

Our previous biochemical studies showed that MLN4924 treatment of AML cell lines and *de novo* AML *ex vivo* results in inhibition of CRL Neddylation, accumulation of the CRL substrate c-Myc, and c-Myc mediated transactivation of the *PMAIP* gene that encodes the proapoptotic Bcl-2 family member Noxa [27, 43]. While performing these studies in *de novo* AML [27], we observed that MLN4924 also induced apoptosis in post-MDS AML (Fig. 3A). Similar to *de novo* AML, post-MDS AML cells treated with MLN4924 exhibited decreased Neddylation of CRLs (Fig. 3B, top panel), increased levels of c-Myc (Fig. 3B, second panel), and upregulation of Noxa at the mRNA and protein levels (Fig. 3B, 3C). Although *de novo* AML and post-MDS AML are distinct diseases that likely occur through different pathologic mechanisms [44], MLN4924 appears to affect both of these forms of AML similarly.

In light of these results, we applied the *ex vivo* drug sensitivity assay to marrow mononuclear cell fractions from two post-MDS AML patients (supplemental online Table 1,


patients 6 and 7). As shown in Figure 2, stem and myeloid progenitor cell populations were reduced in a dose dependent manner in response to increasing concentrations of MLN4924 (Table 1; Fig. 2, patients 6 and 7). Furthermore, populations in these two patients were also less responsive to cytarabine.

To follow up these results, we also assessed the effects of MLN4924 on stem and myeloid progenitor cells from patients with various stages of MDS. Given that these populations contribute to disease initiation and progression before transformation to AML, the ability to target them may also be of therapeutic benefit [1, 9, 45, 46]. Accordingly, we examined samples isolated from five patients with MDS of varying risk category (supplemental online Table 2). Again, the number of events collected, cell count in the starting Live/Lin<sup>-</sup> gate used for analysis (i.e., the population used for analysis, as illustrated in the left panel of Fig. 1A), and the cell count in each stem and progenitor cell population are reported (Table 2).

Previous studies conducted to quantify stem and progenitor cell populations in MDS patients showed that the relative distribution of populations correlates with International Prognostic Scoring System MDS risk category [1, 9, 45, 47]. Two observations comparing low/intermediate risk and high-risk MDS population profiles in our samples were consistent with these reports. First, in the diluent-treated control samples from low and intermediate-risk MDS patients 1 through 4, we noted expansion of the myeloid progenitor cell population relative to stem cell population (Table 2; stem/progenitor cell counts: 6/202, 5/175, 35/190, and 1/98 in samples from patients 1 through 4, respectively). In addition, the MEP populations were suppressed in these low- and intermediate-risk cases (Table 2, patients 1 through 4), another observation previously reported in non-del5q cases in these risk categories [1]. Second, in the sample from high-risk MDS patient 5, we noted a striking expansion of the stem cell compartment relative to the myeloid progenitor cell compartment (Table 2; stem/progenitor cell count: 753/2370). These observations indicated that relative populations observed in our diluent-treated samples after 24 hours of culture *ex vivo* mirrored previous reports regarding alterations in stem and progenitor cell proportions in MDS.

In Figure 4, survival of cells in the MLN4924 treated samples relative to diluent control sample is reported for measurable populations from each patient. MLN4924 almost universally reduces stem and myeloid progenitor cell numbers in MDS, with one exception being the HSCs from the single low-risk MDS patient (Fig. 4B). In this sample, the stem cell counts (HSC + MPP) in the diluent control and MLN4924 treated sample were the same (Table 2, patient 2). Admittedly, these are small populations (five cells). However, reduction in other small populations could readily be seen between diluent and MLN4924 treated samples in other patient samples (Table 2; patient 1, stem cell population and patient 3, GMP population). Moreover, dose-dependent reductions in small populations were noted in the patient 4 samples (Fig. 4D; Table 2).

Although the results are limited by small numbers of stem and progenitor cells in some of the MDS samples, these data collectively demonstrate that drug sensitivity of MDS hematopoietic stem and myeloid progenitor cells can be measured *ex vivo*. Moreover, although the results indicate some variability in sensitivity to MLN4924 among MDS samples, the majority appear to be sensitive to this agent.



**Figure 2.** Survival of AML hematopoietic stem and progenitor cells in MLN4924 or cytarabine versus diluent. Aliquots of bone marrow mononuclear cells from AML patients were treated for 24 hours with increasing concentrations of MLN4924, 100 nM AraC, or diluent (0.1% DMSO). Cells were harvested, stained using the fluorescent antibody panel for identification of live stem and progenitor cells (see Methods), and analyzed by flow cytometry. HSC and progenitor populations were quantified using the gating scheme outlined in Figure 1A. Summary of survival of AML stem and progenitor cells relative to populations in diluent treated samples from seven AML patient BMs. Patient data and cell count in each population are reported in Table 1 and supplemental online Table 1 (de novo AML,  $n = 5$ , patients 1 through 5; post-MDS AML,  $n = 2$ , patients 6 and 7). Several populations could not be gated for quantification from patient 4 nor could the MPP from patient 3 (Table 1).  $CD34^+ CD45^{dim}$  = all stem and myeloid progenitor populations; stem cells = HSC + MPP; myeloid progenitors = CMP + GMP + MEP. Abbreviations: AraC, cytarabine; AML, acute myelogenous leukemia; CMP, common myeloid progenitor; DMSO, dimethyl sulfoxide; GMP, granulocyte-monocyte progenitor; HSC, hematopoietic stem cell; MDS, myelodysplastic syndrome; MEP, megakaryocyte-erythroid progenitor; MPP, multipotent progenitor.


## DISCUSSION

The method we present here for assessing drug sensitivity in MDS and AML stem and progenitor cells contributes a novel dimension to the evaluation of therapeutic agents. Restoring production of functional mature myeloid cells will always be a goal in MDS therapy just as reduction of tumor burden will always be a goal in the treatment of leukemia. However, in the search for agents that produce durable responses, the question is not whether an agent can achieve these immediate goals, but whether it can eliminate the cells responsible for disease initiation, maintenance, and progression. Based on accumulating evidence that stem and progenitor cells play pivotal roles in various stages of MDS and AML pathogenesis [1–3, 7, 9, 48, 49], it has been suggested that targeting one or more of these populations will likely be necessary to eradicate these disorders [5, 6, 13, 16, 18, 19, 50, 51]. This realization provides the rationale for developing a

method to assess the effects of therapeutic agents on MDS and AML stem and progenitor cells.

Several factors influenced the development of the method and studies presented here. First, to diminish cellular stress, decrease processing time, and reduce loss of rare populations that could occur during sorting, the assessment of drug sensitivity in stem and progenitor cells described here (Fig. 1) utilized freshly isolated, unsorted mononuclear cell fractions. Minimizing sample manipulation also reduced the number of variables affecting each sample, allowing changes in stem and progenitor cell number to be attributed to drug exposure. Second, MLN4924 was chosen as a prototype investigational agent based on its continued study in clinical trials for treatment of myeloid disease and our observations of its activity in de novo [27] and post-MDS AML (Fig. 3).

Contrary to normal isolates treated with MLN4924 (Fig. 1C; supplemental online Fig. 1), stem and progenitor cell populations



**Figure 3.** MLN4924 induces apoptosis in post-MDS AML. (**A–C**): Bone marrow mononuclear cells from a patient with post-MDS AML were treated with increasing concentrations of MLN4924 versus diluent (0.1% DMSO) for 24 hours. Parallel samples were assayed for subdiploid (apoptotic) cells by propidium iodide staining (**A**); Cullin1, c-Myc, and Bcl-2 family member expression by immunoblotting (**B**); and Noxa mRNA levels (expressed relative to diluent control) by quantitative RT-PCR (**C**). Note the increase in c-Myc and Noxa as previously seen in de novo AML [27]. Abbreviations: AML, acute myelogenous leukemia; DMSO, dimethyl sulfoxide; GAPDH, glyceraldehyde 3-phosphate dehydrogenase; MDS, myelodysplastic syndrome.

**Table 2.** Cell counts in MDS stem and progenitor populations

Patient	Sample	Total event count	Live/Lin <sup>-</sup> count	34 <sup>+</sup> 45 <sup>dim,a</sup>	Stem cells <sup>b</sup>	HSC	MPP	Myeloid progenitors <sup>c</sup>	CMP	GMP	MEP
1	DMSO	818,238	373,669	209	6	ND	ND	202	ND	ND	ND
	1,000 nM MLN4924	789,594	379,226	124	1	ND	ND	123	ND	ND	ND
2	DMSO	1,677,939	44,702	180	5	2	2	175	98	69	ND
	1,000 nM MLN4924	1,618,407	46,361	111	5	2	2	105	52	45	ND
3	DMSO	1,003,557	21,460	231	35	32	2	190	155	27	6
	1,000 nM MLN4924	1,011,323	22,283	133	23	20	2	106	89	12	0
4	DMSO	1,000,000	1,772	99	1	ND	ND	98	38	43	5
	300 nM MLN4924	1,000,000	2,868	69	ND	ND	ND	69	29	24	2
	1,000 nM MLN4924	1,005,375	3,306	34	ND	ND	ND	34	16	14	0
5	DMSO	727,203	37,606	3,177	753	675	54	2,370	1,974	251	57
	1,000 nM MLN4924	727,608	36,307	2,330	504	504	37	1,788	1,478	183	32

Abbreviations: CD34<sup>+</sup> CD45<sup>dim</sup>, all stem and progenitor populations; CMP, common myeloid progenitor; DMSO, dimethyl sulfoxide; GMP, granulocyte-monocyte progenitor; HSC, hematopoietic stem cell; Lin<sup>-</sup>, Lineage negative; MDS, myelodysplastic syndrome; MEP, megakaryocyte-erythroid progenitor; MPP, multipotent progenitor; ND, population not detected.


<sup>a</sup>CD34<sup>+</sup> CD45<sup>dim</sup> = starting population for analysis of all subpopulations (Fig. 1A).

<sup>b</sup>Stem cells = HSC + MPP.

<sup>c</sup>Myeloid progenitors = CMP + GMP + MEP.

in de novo and secondary AML were reduced by MLN4924 treatment ex vivo (Fig. 2). Exceptions to this trend occurred in HSC, MPP, CMP, and MEP populations from FLT3<sup>+</sup> AML patient 1 and the HSC and MPP populations from patient 2

(Fig. 2). Although previous studies in our laboratory and others reported cytotoxic effects of MLN4924 in FLT3 mutated AML cell lines and bulk leukemic cell isolates from patients [27, 52], results from the present study suggest that



**Figure 4.** Quantification and relative survival of MDS hematopoietic stem and progenitor cells treated with MLN4924 versus diluent. Bone marrow mononuclear cells were treated for 24 hours with MLN4924 or diluent (0.1% DMSO), harvested, stained using the fluorescent antibody panel for identification of live stem and progenitor cells (see Methods) and analyzed by flow cytometry. HSC and myeloid progenitor populations were quantified using the gating scheme outlined in Figure 1A. (A–E): Survival for each population in samples treated with MLN4924 relative to corresponding diluent-treated population for MDS patients 1 through 5. Patient data and cell count in each population are reported in Table 2 (Figure legend continues on next page.)



FLT3<sup>+</sup> AML stem and progenitor cells may not be as sensitive to MLN4924. Given the need for additional therapeutic options in this patient group, future study of additional samples is warranted. In contrast to this discrepancy between bulk cell population and stem and progenitor cell drug sensitivity assays, the cytotoxic effect of MLN4924 seen in the post-MDS AML bulk cell isolate (Fig. 3A) was also observed in post-MDS AML stem and progenitor cells from AML patients 6 and 7 (closed square and open circle, respectively, Fig. 2). Also of note in these two samples, cell counts in HSC, MPP, CMP, GMP, and MEP populations were all reduced in MLN4924 treated aliquots compared with diluent control. Although only two samples were assessed, this result was especially intriguing given that post-MDS AML patients do not respond optimally to standard cytarabine-based chemotherapy, with median overall survival of only 6.5 months [53].

Another trend observed in de novo AML stem and progenitor cells is the lack of response to cytarabine (Fig. 2). Although detailed investigation of the mechanisms underlying the observed differences between MLN4924 and cytarabine is beyond the scope of this study, it is clear that the leukemic stem and progenitor cell response to these agents differs. The ex vivo cytarabine exposure utilized in the present study is considerably shorter than the 7-day treatment often used during induction in vivo and the concentration used is lower than the cytarabine levels transiently achieved during consolidation therapy [54]. Therefore, it is possible that the response of leukemic stem and progenitor cells to cytarabine is greater in vivo. Further investigation is required to assess this possibility.

Cells isolated from low- and high-risk MDS patients (supplemental online Table 2; Table 2) exhibited the relative distributions of stem and progenitor cells previously reported in these two MDS patient groups [1, 9, 45]. Specifically, in diluent-treated samples from low- and intermediate-risk MDS (patients 1 through 4), greater progenitor cell counts relative to stem cell counts were observed, whereas in high-risk MDS (patient 5), the stem cell population was expanded relative to the progenitor cell population (supplemental online Table 2; Table 2). These observations provide evidence that the MDS samples examined in the present study are representative of this disorder.

In studies performed on parallel aliquots from these MDS patient samples (supplemental online Table 2), MLN4924 treatment reduced the number of MDS stem (HSCs and MPPs) and progenitor cells (CMPs, GMP, and MEP) compared with their respective numbers present in the diluent control in all five MDS patient cell isolates (Fig. 4). One exception was noted in the stem cell population (HSC + MPP) from patient 2 (Fig. 4B). To our knowledge, this is the first report of MLN4924 activity in MDS samples before leukemic transformation. These results by themselves cannot conclusively establish that only MDS stem and progenitor cells, as opposed to “normal” stem cells, were affected. However, even in low risk MDS patients,

cytogenetic and mutation profiling studies have previously demonstrated that the overwhelming majority of HSCs and progenitors are abnormal [45], making it probable that MLN4924 reduction of these populations includes the mutated HSCs and progenitors. Moreover, our study as well as a recently published study from another group [29] indicate that normal stem and progenitor cells are unaffected by the same concentrations of MLN4924 (Fig. 1C; supplemental online Fig. 1).

Although MDS and AML are notoriously heterogeneous diseases, the contribution of malignant stem and progenitor cells to disease initiation and progression is a unifying feature [1, 2, 44, 48, 55]. Defining which specific stem or progenitor cell populations need to be targeted to induce durable remission is an active area of investigation [18]. Furthermore, response or resistance of these malignant stem and progenitor populations has been shown to correlate with clinical course and outcome [1, 2, 6, 8, 56]. Although the strategy presented here was used to study the effects of cytarabine and MLN4924, this method provides an individualized, patient specific approach for ex vivo drug sensitivity testing in leukemic stem and progenitor cell populations that could be applied to other agents. In combination with other tumor profiling platforms and clinical parameters, this approach could potentially be utilized to inform design and selection of therapeutic agents.

## CONCLUSION

The method presented here enables assessment of drug response in malignant hematopoietic stem and progenitor cells. Because these cells play a critical role in disease initiation, maintenance, and progression, the response of these cells is of special interest when considering the ability of an agent to produce durable clinical responses.

## ACKNOWLEDGMENTS

This work was supported in part by NIH Grants F30 CA183507 (K.L.B.K.), T32 GM65841 (K.L.B.K.), and R01 CA166704 (S.H.K.) as well as P30 CA06973 and U01 CA70095 (to support sample acquisition). Additional support was provided through the Carol Bacon Bone Marrow Transplant and Leukemia Research Award (L.E.F. and K.L.B.K.). We thank the expert staff members of the Mayo Clinic Rochester Flow Cytometry Core for their advice and assistance, Jacqueline Joy and Jackie Greer for coordination of sample acquisition, and Deb Strauss for administrative support.

## AUTHOR CONTRIBUTIONS:

K.L.B.K.: conception and design, financial support, collection and/or assembly of data, data analysis and interpretation, manuscript writing, final approval of manuscript; L.E.F.: conception and design, financial support, provision of study material or patients,

(Figure legend continued from previous page.)

and supplemental online Table 2. CD34<sup>+</sup> CD45<sup>dim</sup> = all stem and myeloid progenitor populations; stem cells = HSC + MPP; myeloid progenitors = CMP + GMP + MEP. Abbreviations: CMP, common myeloid progenitor; DMSO, dimethyl sulfoxide; GMP, granulocyte-monocyte progenitor; HSC, hematopoietic stem cell; MEP, megakaryocyte-erythroid progenitor; MPP, multipotent progenitor; ND, population not detected.

data analysis and interpretation, final approval of manuscript; B.D.S., A.D.H., and J.M.F.: provision of study material or patients, final approval of manuscript; J.E.K.: provision of study material or patients, data analysis and interpretation, final approval of manuscript; S.H.K.: conception and design, financial support, data analysis and interpretation, manuscript writing, final approval of manuscript.

#### DISCLOSURE OF POTENTIAL CONFLICTS OF INTEREST

J.M.F. has research support for the C15009 clinical trial (azacitidine & MLN4924, Millennium/Takeda). J.E.K. is a consultant for Tolero Pharmaceuticals. S.H.K. has research funding from Eli Lilly for an unrelated project. The other authors indicated no potential conflicts of interest.

#### REFERENCES

- Woll PS, Kjällquist U, Chowdhury O et al. Myelodysplastic syndromes are propagated by rare and distinct human cancer stem cells in vivo. *Cancer Cell* 2014;25:794–808.
- Shlush LI, Zandi S, Mitchell A et al. Identification of pre-leukaemic haematopoietic stem cells in acute leukaemia. *Nature* 2014;506:328–333.
- Lapidot T, Sirard C, Vormoor J et al. A cell initiating human acute myeloid leukaemia after transplantation into SCID mice. *Nature* 1994;367:645–648.
- Pollyea DA, Gutman JA, Gore L et al. Targeting acute myeloid leukemia stem cells: A review and principles for the development of clinical trials. *Haematologica* 2014;99:1277–1284.
- Kreso A, Dick JE. Evolution of the cancer stem cell model. *Cell Stem Cell* 2014;14:275–291.
- Craddock C, Quek L, Goardon N et al. Azacitidine fails to eradicate leukemic stem/progenitor cell populations in patients with acute myeloid leukemia and myelodysplasia. *Leukemia* 2013;27:1028–1036.
- Welch JS, Ley TJ, Link DC et al. The origin and evolution of mutations in acute myeloid leukemia. *Cell* 2012;150:264–278.
- Goardon N, Marchi E, Atzberger A et al. Coexistence of LMPP-like and GMP-like leukemia stem cells in acute myeloid leukemia. *Cancer Cell* 2011;19:138–152.
- Pang WW, Pluvinau JV, Price EA et al. Hematopoietic stem cell and progenitor cell mechanisms in myelodysplastic syndromes. *Proc Natl Acad Sci USA* 2013;110:3011–3016.
- Taussig DC, Miraki-Moud F, Anjos-Afonso F et al. Anti-CD38 antibody-mediated clearance of human repopulating cells masks the heterogeneity of leukemia-initiating cells. *Blood* 2008;112:568–575.
- Tomasetti C, Vogelstein B. Cancer etiology. Variation in cancer risk among tissues can be explained by the number of stem cell divisions. *Science* 2015;347:78–81.
- Liesveld J. Management of AML: Who do we really cure? *Leuk Res* 2012;36:1475–1480.
- Krause DS, Van Etten RA. Right on target: Eradicating leukemic stem cells. *Trends Mol Med* 2007;13:470–481.
- Lichtman MA. Myelodysplasia or myeloneoplasia: Thoughts on the nosology of clonal myeloid diseases. *Blood Cells Mol Dis* 2000;26:572–581.
- Lichtman MA. The stem cell in the pathogenesis and treatment of myelogenous leukemia: A perspective. *Leukemia* 2001;15:1489–1494.
- Dick JE. Stem cell concepts renew cancer research. *Blood* 2008;112:4793–4807.
- Clarke MF, Dick JE, Dirks PB et al. Cancer stem cells—perspectives on current status and future directions: AACR Workshop on cancer stem cells. *Cancer Res* 2006;66:9339–9344.
- Meacham CE, Morrison SJ. Tumour heterogeneity and cancer cell plasticity. *Nature* 2013;501:328–337.
- Hanahan D, Weinberg RA. Hallmarks of cancer: The next generation. *Cell* 2011;144:646–674.
- Kufe DW, Major PP, Egan EM et al. Correlation of cytotoxicity with incorporation of ara-C into DNA. *J Biol Chem* 1980;255:8997–9000.
- Kufmann SH. Induction of endonucleolytic DNA cleavage in human acute myelogenous leukemia cells by etoposide, camptothecin, and other cytotoxic anticancer drugs: A cautionary note. *Cancer Res* 1989;49:5870–5878.
- Appelbaum FR, Rowe JM, Radich J et al. Acute myeloid leukemia. *Hematology (Am Soc Hematol Educ Program)* 2001;2001:62–86.
- Tallman MS, Gilliland DG, Rowe JM. Drug therapy for acute myeloid leukemia. *Blood* 2005;106:1154–1163.
- Brownell JE, Sintchak MD, Gavin JM et al. Substrate-assisted inhibition of ubiquitin-like protein-activating enzymes: The NEDD8 E1 inhibitor MLN4924 forms a NEDD8-AMP mimetic in situ. *Mol Cell* 2010;37:102–111.
- Soucy TA, Smith PG, Milhollen MA et al. An inhibitor of NEDD8-activating enzyme as a new approach to treat cancer. *Nature* 2009;458:732–736.
- Mathewson N, Toubai T, Kapeles S et al. Neddylation plays an important role in the regulation of murine and human dendritic cell function. *Blood* 2013;122:2062–2073.
- Knorr KL, Schneider PA, Meng XW et al. MLN4924 induces Noxa upregulation in acute myelogenous leukemia and synergizes with Bcl-2 inhibitors. *Cell Death Differ* 2015;22:2133–2142.
- Swords RT, Erba HP, DeAngelo DJ et al. Pevonedistat (MLN4924), a first-in-class NEDD8-activating enzyme inhibitor, in patients with acute myeloid leukaemia and myelodysplastic syndromes: A phase 1 study. *Br J Haematol* 2015;169:534–543.
- Zhou L, Chen S, Zhang Y et al. The NAE inhibitor pevonedistat interacts with the HDAC inhibitor belinostat to target AML cells by disrupting the DDR. *Blood* 2016;127:2219–2230.
- English D, Andersen BR. Single-step separation of red blood cells. Granulocytes and mononuclear leukocytes on discontinuous density gradients of Ficoll-Hypaque. *J Immunol Methods* 1974;5:249–252.
- Rocha V, Labopin M, Sanz G et al. Transplants of umbilical-cord blood or bone marrow from unrelated donors in adults with acute leukemia. *N Engl J Med* 2004;351:2276–2285.
- Gluckman E, Broxmeyer HA, Auerbach AD et al. Hematopoietic reconstitution in a patient with Fanconi's anemia by means of umbilical-cord blood from an HLA-identical sibling. *N Engl J Med* 1989;321:1174–1178.
- Bender JG, Unverzagt K, Walker DE et al. Phenotypic analysis and characterization of CD34+ cells from normal human bone marrow, cord blood, peripheral blood, and mobilized peripheral blood from patients undergoing autologous stem cell transplantation. *Clin Immunol Immunopathol* 1994;70:10–18.
- Van Epps DE, Bender J, Lee W et al. Harvesting, characterization, and culture of CD34+ cells from human bone marrow, peripheral blood, and cord blood. *Blood Cells* 1994;20:411–423.
- Basford C, Forraz N, McGuckin C. Optimized multiparametric immunophenotyping of umbilical cord blood cells by flow cytometry. *Nat Protoc* 2010;5:1337–1346.
- Nicoletti I, Migliorati G, Pagliacci MC et al. A rapid and simple method for measuring thymocyte apoptosis by propidium iodide staining and flow cytometry. *J Immunol Methods* 1991;139:271–279.
- Kaufmann SH, Lee SH, Meng XW et al. Apoptosis-associated caspase activation assays. *Methods* 2008;44:262–272.
- Yun S, Vincelette ND, Knorr KLB et al. 4EBP1/c-MYC/PUMA and NFB/EGR1/BIM pathways underlie cytotoxicity of mTOR dual inhibitors in malignant lymphoid cells. *Blood* 2016;127:2711–2722.
- Caserta TM, Smith AN, Gultice AD et al. Q-VD-OPh, a broad spectrum caspase inhibitor with potent antiapoptotic properties. *Apoptosis* 2003;8:345–352.
- Milhollen MA, Traore T, Adams-Duffy J et al. MLN4924, a NEDD8-activating enzyme inhibitor, is active in diffuse large B-cell lymphoma models: Rationale for treatment of NF-kappaB-dependent lymphoma. *Blood* 2010;116:1515–1523.
- Walters RS, Kantarjian HM, Keating MJ et al. Mitoxantrone and high-dose cytosine arabinoside in refractory acute myelogenous leukemia. *Cancer* 1988;62:677–682.
- Heinemann V, Murray D, Walters R et al. Mitoxantrone-induced DNA damage in leukemia cells is enhanced by treatment with high-dose arabinosylcytosine. *Cancer Chemother Pharmacol* 1988;22:205–210.
- Nikiforov MA, Riblett M, Tang WH et al. Tumor cell-selective regulation of NOXA by c-MYC in response to proteasome inhibition. *Proc Natl Acad Sci USA* 2007;104:19488–19493.
- Cazzola M, Della Porta MG, Malcovati L. The genetic basis of myelodysplasia and its clinical relevance. *Blood* 2013;122:4021–4034.
- Will B, Zhou L, Vogler TO et al. Stem and progenitor cells in myelodysplastic syndromes show aberrant stage-specific expansion and harbor genetic and epigenetic alterations. *Blood* 2012;120:2076–2086.
- Porwit A, van de Loosdrecht AA, Bettelheim P et al. Revisiting guidelines for

integration of flow cytometry results in the WHO classification of myelodysplastic syndromes-proposal from the International/European LeukemiaNet Working Group for Flow Cytometry in MDS. *Leukemia* 2014;28:1793–1798.

**47** van de Loosdrecht AA, Ireland R, Kern W et al. Rationale for the clinical application of flow cytometry in patients with myelodysplastic syndromes: Position paper of an International Consortium and the European LeukemiaNet Working Group. *Leuk Lymphoma* 2013;54:472–475.

**48** Jan M, Snyder TM, Corces-Zimmerman MR et al. Clonal evolution of preleukemic hematopoietic stem cells precedes human acute myeloid leukemia. *Sci Transl Med* 2012;4:149ra118.

**49** Corces-Zimmerman MR, Hong WJ, Weissman IL et al. Preleukemic mutations in human acute myeloid leukemia affect epigenetic regulators and persist in remission. *Proc Natl Acad Sci USA* 2014;111:2548–2553.

**50** Chen WC, Yuan JS, Xing Y et al. An integrated analysis of heterogeneous drug responses in acute myeloid leukemia that enables the discovery of predictive biomarkers. *Cancer Res* 2016;76:1214–1224.

**51** Tehranchi R, Woll PS, Anderson K et al. Persistent malignant stem cells in del(5q) myelodysplasia in remission. *N Engl J Med* 2010;363:1025–1037.

**52** Swords RT, Kelly KR, Smith PG et al. Inhibition of NEDD8-activating enzyme: A novel approach for the treatment of acute myeloid leukemia. *Blood* 2010;115:3796–3800.

**53** Bello C, Yu D, Komrokji RS et al. Outcomes after induction chemotherapy in patients with acute myeloid leukemia arising from myelodysplastic syndrome. *Cancer* 2011;117:1463–1469.

**54** Donehower RC, Karp JE, Burke PJ. Pharmacology and toxicity of high-dose cytarabine by 72-hour continuous infusion. *Cancer Treat Rep* 1986;70:1059–1065.

**55** Jaiswal S, Fontanillas P, Flannick J et al. Age-related clonal hematopoiesis associated with adverse outcomes. *N Engl J Med* 2014;371:2488–2498.

**56** Griessinger E, Anjos-Afonso F, Vargaftig J et al. Frequency and dynamics of leukemia-initiating cells during short-term ex vivo culture informs outcomes in acute myeloid leukemia patients. *Cancer Res* 2016;76:2082–2086.



See [www.StemCellsTM.com](http://www.StemCellsTM.com) for supporting information available online.